
Dr. Carlo Marassi, DDS, MSc.
University of São Paulo - Bauru, SLM Research Center
Private Practice - Rio de Janeiro - Brazil

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1. Cuperus AM, Harms MC, Rangel FA, Bronkhorst EM, Schols JG, Breuning KH. Dental models made with an intraoral scanner: a validation study. *Am J Orthod Dentofacial Orthop.* 2012;142(3):308-13.
2. Park JH, Laslovich J. Trends in the Use of Digital Study Models and Other Technologies Among Practicing Orthodontists. *J Clin Orthod.* 2016;50(7):413-9.
3. Burzynski JA, Firestone AR, Beck FM, Fields HW, Deguchi T. Comparison of digital intraoral scanners and alginate impressions: Time and patient satisfaction. *Am J Orthod Dentofacial Orthop.* 2018;153(4):534-41.
4. Fleming PS, Marinho V, Johal A. Orthodontic measurements on digital study models compared with plaster models: a systematic review. *Orthod Craniofac Res.* 2011;14(1):1-16.
5. Breuning KH, Kau CH. *Digital Planning and Custom Orthodontic Treatment.* 2017
6. Lombardo, L.; Arreghini, A.; Ramina, F.; Huanca Ghislanzoni, L.; Siciliani, G. Predictability of orthodontic movement with orthodontic aligners: A retrospective study. *Prog. Orthod.* 2017, 18, 35.
7. Zheng M, Liu R, Ni Z, Yu Z. Efficiency, effectiveness and treatment stability of clear aligners: A systematic review and meta-analysis. *Orthod Craniofac Res.* 2017;20:127–33.
8. Ke Y, Zhu Y, Zhu M, A comparison of treatment effectiveness between clear aligner and fixed appliance therapies. *BMC oral health* 2019; 19(1), 24.
9. Rossini G, Parrini S, Castrolforio T, Deregibus A, Debernardi CL. Efficacy of Clear Aligners in Controlling Orthodontic Tooth Movement : A Systematic Review. *Angle Orthod.* 2015;85:881–9

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About Dr. Carlo Marassi

Graduated cum laude at Federal University - Rio de Janeiro Specialization in Orthodontics at State University - São Paulo Master of science at SLM Research Center - Campinas Emeritus Editor - Journal of Clinical Orthodontics - Dental Press Former Scientific director of the Brazilian Society of Orthodontists Vice director of the Society of Orthodontics of Rio de Janeiro Scientific director of the Straight-Wire Group of Rio de Janeiro Residence - Michigan University - Ann Arbor - United States President and CEO of Marassi Training Institute International Lecturer

About Dr. Alice Spitz

Graduated at Federal University of Rio de Janeiro Specialization in Orthodontics at Federal University of Rio de Janeiro Master of science at Federal University of Rio de Janeiro PHD at Federal University of RJ / University of Pennsylvania - USA Collaborating Professor - Department of Orthodontics - Federal University of Bahia

About 3Shape

3Shape is changing dentistry together with dental professionals across the world by developing innovations that provide superior dental care for patients. Our portfolio of 3D scanners and CAD/CAM software solutions for the dental industry includes the multiple award-winning 3Shape TRIOS® intraoral scanner, the upcoming 3Shape X1® CBCT scanner, as well as market-leading scanning and design software solutions for both dental practices and labs.

Two graduate students founded 3Shape in Denmark's capital in the year 2000. Today, 3Shape has over 1,500 employees serving customers in over 100 countries from 3Shape offices around the world. 3Shape's products and innovations continue to challenge traditional methods, enabling dental professionals to treat more patients more effectively. www.3shape.com

Let's change dentistry together

www.3shape.com

3shape 