

SPARK™

CLEAR ALIGNER SYSTEM

3SHAPE INTEGRATION QUICK START GUIDE



SPARK QUICK START GUIDE


SETTING UP A CONNECTION WITH 3SHAPE

Spark doctors who use 3Shape scanners now can more quickly and easily access and submit scans for Spark patients through the Spark Aligner website. This document reviews the simple steps to set up a one time connection with their 3 Shape Communicate™ account through the Spark Aligner website saving time and effort.

Please ensure you have your 3Shape Communicate Portal email and password before using this guide.

STEP 1

Login to the Spark Aligner Website at aligners.ormcodigital.com.



STEP 2

Under your doctor profile menu, select **My Profile**.

The screenshot shows the Ormco | DTX user interface. At the top right, there is a user profile dropdown menu for 'Sample Doctor'. The menu is open, showing three options: 'My Profile', 'Clinical Preferences', and 'Logout'. A red arrow points to the 'My Profile' option. Below the menu, there is a search bar and a table of patient cases.

Patient	Status	Submitted	Order	Office Notes
104049 Patient A	Unsubmitted Submit Case		Primary	Add Note
105512 Patient B	Unsubmitted Submit Case		Primary	Add Note
103028 Patient C	Awaiting Approval Open Approver		Primary	Add Note
109503 Patient D	On Hold Contact Support		Primary Spark 20	Add Note

STEP 3

Scroll down to the **Impression Submission Preferences** section and select **3Shape Scans**.

The screenshot shows the 'My Profile - Doctor' page in the Ormco | DTX system. The 'Impression Submission Preferences' section is visible. Under 'Impression Type:', there are three radio button options: 'Scans / STL', 'iTero Scans', and '3Shape Scans'. The '3Shape Scans' option is selected, and a red arrow points to it. Below the selected option, the status is 'Not Connected' and there is a link to 'Connect 3Shape to Spark'.

STEP 4

Click on the **Connect 3Shape to Spark** link.

The screenshot shows the Ormco | DTX user interface. At the top, there is a navigation bar with 'Ormco | DTX' on the left and 'Patients', 'Help', and 'Sample Doctor' on the right. Below the navigation bar, there is a 'My Profile- Doctor' section with a note: 'Note: Fields marked with an asterisk* are required'. Underneath is the 'Impression Submission Preferences' section. It contains 'Impression Type' with three radio button options: 'Scans / STL', 'iTero Scans', and '3Shape Scans' (which is selected). Below the radio buttons, it says 'Status: Not Connected' and a blue link labeled 'Connect 3Shape to Spark'. A red arrow points to this link.


STEP 5

Upon clicking, you will be taken directly to the 3Shape Communicate Portal. Please enter your email/username and password. Then click **Sign In**.

The screenshot shows a web browser window with the title 'Log in - 3Shape Identity'. The browser address bar shows a secure connection. The page features the 3shape logo in the top left. The main heading is 'Log in' with the subtext 'With your 3Shape account'. There are two input fields: 'EMAIL/USERNAME' containing 'testOrmcoLabUS@3shapecommunicate.com' and 'PASSWORD' with masked characters. Below the password field is a red link for 'Forgot password?'. At the bottom, there is a prominent red 'Sign in' button.


STEP 6

Upon signing in, you will be automatically redirected back to the Spark Aligner Website when your connection has been made. Click **OK** to be taken back to your profile page on the Spark Aligner Website.



STEP 7

Scroll down your profile page to ensure that it is showing **Connected**.



IMPORTANT NOTE

Doctors still need to ensure that Ormco lab connection is made within their 3Shape Communicate Portal and that all scans for Spark patients are sent to Ormco during the scanning process.

Go to <https://portal.3shapecommunicate.com/connections>

Click on **Add Connection**




Search for the relevant Ormco lab account based on your region:

EU: Insignia.SupportEurope@ormco.com

US: Insignia.Support@ormco.com

ASIA: Insignia.SupportAsia@ormco.com

Click **Connect**





NEED HELP?

Please contact Spark Customer Support at:

North America

Email: spark.support@ormco.com

Call: 714.516.7400

800.854.1741

Australia & New Zealand

Email: ormco.australia@sybrondental.com

Call: 1800 023 603

Europe

Email: spark.csdigital@ormco.com

Call: 00800 3032 3032

900 983 191