

For Immediate Release

Contact: Bruce Frederic Mendel, 3Shape Communications

E-mail: bruce.mendel@3shape.com

Web: <u>www.3shape.com</u>

Date: September 6, 2023 **HQ Copenhagen, Denmark**

Engaging patients with digital dentistry has never been easier

- New and improved award-winning MOVE+ cart released
- Compatible with all TRIOS intraoral scanner models
- Can be purchased with or without a scanner.

Copenhagen, September 6, 2023 – 3Shape today announced the release of a new and improved version of its award-winning MOVE+ cart for TRIOS intraoral scanners.

The new MOVE+ cart makes it easier than ever for dentists to engage and build trust with their patients.

Dental professionals can take advantage of the MOVE+ award-winning ergonomic design to easily share digital scans, images, and treatment plans with patients. Because the MOVE+ cart is on wheels and its swivel arm holding the 15.6" touch screen is fully adjustable, professionals can glide the scanner, screen, and stand wherever they want in the operatory.

MOVE+ easy ergonomics ensure that the screen and scanner are always in the right place to engage, excite, and drive treatment acceptance with patients.

"We know patients are excited and impressed when they see their intraoral scans. And we know that dental professionals love sharing scans with their patients because it makes explaining treatment more transparent and impactful. The new MOVE+ makes it so much easier to do both," says Rune Fisker, 3Shape Senior Vice President for Product Strategy.

He adds, "Dentists are always looking for creative ways to position the PC screen and move their intraoral scanner around the operatory. MOVE+ solves these challenges and does so in a Red Dot award-winning elegant cart design."

The new and improved MOVE+ has a more powerful Windows 11 PC, making scan post-processing time 30% faster. Its premium brushed aluminum esthetic contours complement any dental practice's interior design.

All 3Shape TRIOS models are compatible with MOVE+, including the recently released TRIOS 5. MOVE+ is available with or without a TRIOS intraoral scanner. The stand connects to a practice's internet via Wi-Fi or fixed cable.



For Immediate Release

MOVE+ specs:

Intel® Core™ i5-12600T Windows® 11 IoT Enterprise Value, 64bit 1TB SSD storage T1000 8GB GDDR6 GPU 32GB DDR4 RAM Weighs 14kg (30.8 lbs.) MSR: 6,000 €

About 3Shape

3Shape was founded in 2000 by two people determined to unlock opportunities in 3D digital technology. From two, we've grown to over 2,000 – and from our beginning in Denmark, we're now present in over 100 countries around the world.

Our singular purpose remains the same – to advance and connect medical professionals and patients. We drive progress in digital dentistry with groundbreaking, industry-recognized scanners for clinics and labs. And connect dentists and technicians with software workflows that just flow. We advance and coach professionals on every step of their digital journey – helping to truly elevate the customer experience.

Our impact is reflected in the results we achieve. Every second of every day, our award-winning products restore the confidence of a patient. Decisively going beyond products, we bring ongoing training and education that not only helps professionals stay relevant, it puts them ahead of the game.

We actively advance and connect dentists and lab technicians – the people who are the backbone of the industry. Advance. Connect. Together.

www.3shape.com





For Immediate Release





For Immediate Release

