

Contact: Bruce Frederic Mendel, 3Shape Communications
E-mail: bruce.mendel@3shape.com
Web: www.3shape.com
Date: June 6, 2016, HQ Copenhagen, Denmark

New Dental System 2016 to Boost Dental Lab Productivity

Copenhagen, June 6, 2016 – 3Shape, a global leader in 3D scanners and dental CAD/CAM software solutions, releases Dental System™ 2016 – the latest version of its industry leading CAD software for dental professionals.

Dental System™ 2016 introduces new and improved applications and workflows as well as a completely rewritten software for impression scanning to eliminate steps in the dental lab workflow. The new Dental System 2016 also saves time for lab technicians by making designing faster and more intuitive for simpler to more advanced cases as well as increasing software reliability.

“Working closely with industry leaders from around the world has helped 3Shape to simplify and create a more efficient, cost effective, and predictable workflow for dental lab professionals,” says Flemming Thorup, president & CEO at 3Shape. “Their expert feedback has been essential in developing an even more powerful and intuitive 2016 version of our Dental System CAD solution.”

Dental System 2016 introduces new features like ground-breaking impression scanning. The 2016 version of the impression scanning application eliminates a dental lab’s need for pouring gypsum models and other manual steps in the workflow.

All 3Shape dental lab scanners can now scan analog impressions and skip directly to CAD design: a workflow step-saver that can reduce crown turnaround time by nearly one day. The 3Shape D2000 and D1000 lab scanner models however, feature multi-line scanning technology - for a more-detailed and optimal impression scanning.

"3Shape’s impression scanning is today's reality and the future. No waiting for models to set, no plaster inaccuracies and we reduce costs. Every lab should invest in this technology," says Carlo Backes, at Excent Tandtechnik BV in the Netherlands.

3Shape Dental System™ 2016 highlights:

- Ground-breaking impression scanning software - Scan each side of a triple tray in less than 30 seconds and single trays in 75 seconds.
- New functionality accelerates 3Shape D500, D1000 and D2000 lab scanners performance and speed
- Increased stability for large cases such as implant bridges

- Integration with TRIOS® and 3rd party IOS
- Robust crown design engine for fast and intuitive design
- RealView™ with 2D Overlay
- Easily work on advanced cases such as RPDs and full dentures

Upgrading to Dental System 2016 is a part of 3Shape LABcare™

3Shape includes LABcare™ with all Dental System subscriptions. LABcare ensures 3Shape users that they continuously receive the latest technology through annual software releases such as Dental System 2016. In addition to upgrades, LABcare includes access to the 3Shape support network, as well as training channels such as webinars, videos and more.

Dental System training video series on YouTube:

https://www.youtube.com/user/3shape/playlists?sort=dd&shelf_id=4&view=50

How to install and upgrade to Dental System 2016:

- How to install Dental System 2016 - https://www.youtube.com/watch?v=_gkrlYJaXUU
- How to upgrade to Dental System 2016 - <https://www.youtube.com/watch?v=DJmyatnJeJU>

Dental System 2016 is available through 3Shape resellers. Availability to end-users is dependent on the specific system configuration. Please contact your local 3Shape supplier, or visit <http://www.3shape.com> regarding reseller information. For more information regarding Dental System 2016: <http://www.3shape.com/new+products/dental+labs/cad+software>

About 3Shape

3Shape creates 3D scanning and CAD/CAM software solutions. Award-winning technology that enables dental and hearing professionals to treat more people, more effectively and with improved care. A privately-owned company, 3Shape has over 700 employees with a product-development force of more than 275 professionals. Offices and service centers located in the Americas, Asia and Europe serve customers in more than 100 countries. Company headquarters are in Copenhagen, Denmark. www.3shape.com.



3Shape Headquarters
Europe, Middle East & Africa
Holmens Kanal 7
1060 Copenhagen, Denmark
Tel: +45 7027 2620

3Shape Asia
Room 906, Tower A of Eton Place
No. 69, Dongfang Road
200120 Shanghai, China
Tel: +86 21 5835 2281

3Shape Latin America
Carrera 13 # 82-91
Oficina 401
110221 Bogotá, Colombia
Tel: +57 1691 9508

3Shape North America
Somerset Hills Corporate Center
10 Independence Boulevard, Suite 150
Warren, New Jersey 07059, USA
Tel: +1 (908) 867 0144

More information at
www.3shape.com
info@3shape.com