PRESS RELEASE



For Immediate Release

Contact: Bruce Frederic Mendel, 3Shape Communications

E-mail: bruce.mendel@3shape.com

Web: <u>www.3shape.com</u>

Date: September 1, 2020 **HQ Copenhagen, Denmark**

Stay up to speed with the brand-new 20% faster Generation Red E scanners from 3Shape

Copenhagen, September 1, 2020 – 3Shape announces the launch of its brand-new Generation Red E scanners for labs.

The Generation Red E scanners are upgrades of the 3Shape E1, E2 and E3 dental lab scanners to the same next generation platform as the award-winning E4 lab scanner. The next generation platform enables the three Generation Red E scanners to scan 20% faster for optimized productivity.

The new Generation Red E scanners also add a signature red ring to the Red Dot award-winning E scanner design.

20% faster at no extra cost!

Despite the significant boost in scanning speed, all Generation Red E scanners are available at the same price points as their corresponding E models.

"With the 20% increase in speed, coupled with the already lightning fast E4, the Generation Red line of lab scanners make it easier than ever for labs to do more of what they love, faster. Productivity is paramount for labs. The speed and accuracy of Generation Red E scanners combined with our easy to use software ensure that lab professionals can work quickly, predictably, and at a great price!" says Rune Fisker, 3Shape Senior Vice President, Product Strategy.

SAVE up to USD \$2,400/€2,000 on a brand-new Generation Red Scanner today

Up until the end November 2020, labs can save up to \$1,200/€1,000 on a new Generation Red E Scanner – and up to \$2,400/€2,000 by trading up their old 3Shape lab scanner for a new Generation Red E scanner. Professionals are encouraged to contact their 3Shape reseller for complete details on the offer.

Generation Red E scanner speeds and accuracy

Generation Red E scanners include from two to four 5MP cameras (depending on the model) and auto start scanning. ISO-documented accuracy in the E4 and E3 models is suitable for implant bars. The E line of lab scanners are Red Dot international product design winners.

	Generation Red	Generation Red	Generation Red	Generation Red
	E4	E3	E2	E1
Full arch scan speed	9 sec.	18 sec.	24 sec.	32 sec.

PRESS RELEASE



For Immediate Release

Full arch impression	45 sec.	64 sec.	72 sec.	104 sec.
scan speed				
Accuracy	4 μm	7 μm	10 μm	10 μm
(ISO 12836)				!

Discover more about the Generation Red E scanners: https://www.3shape.com/en/scanners/lab

Compatibility with the Generation Red E scanners and 3Shape orthodontic software coming in October 2020.

Contact your reseller regarding availability of 3Shape products in your region.

About 3Shape

3Shape is changing dentistry together with dental professionals across the world by developing innovations that provide superior dental care for patients. Our portfolio of 3D scanners and CAD/CAM software solutions includes the multiple award-winning 3Shape TRIOS intraoral scanner, the upcoming 3Shape X1 CBCT scanner, as well as market-leading scanning and design software solutions for both dental practices and labs.

Two graduate students founded 3Shape in Denmark's capital in the year 2000. Today, 3Shape has over 1,600 employees serving customers in over 100 countries from an ever-growing number of 3Shape offices around the world. 3Shape's products and innovations continue to challenge traditional methods, enabling dental professionals to treat more patients more effectively. www.3shape.com



For Immediate Release

