

Contact: Bruce Frederic Mendel, 3Shape Communications
E-mail: bruce.mendel@3shape.com
Web: www.3shape.com
Date: February 23, 2022 HQ Copenhagen, Denmark

From 3Shape TRIOS intraoral scans to dental model designs in minutes

3 new ways 3Shape AI & cloud solutions make digital dentistry more efficient

Copenhagen, February 23, 2022 – 3Shape to introduce its brand-new Model Maker app, Dental Practice Management Software integrations, and 3Shape Automate innovation at LAB DAY 2022.

3Shape Senior Vice President of Product Strategy, Rune Fisker, will present the three new solutions and more during his LAB DAY opening press conference on February 24 at 9.30AM CST – booth J-16 K-17.

What is Model Maker?

Model Maker enables both 3Shape TRIOS intraoral scanner and 3Shape Studio apps users to instantly convert digital impressions into dental models. Model Maker's automated workflow is powered by AI and cloud technology. And converts TRIOS scans to models in just minutes. The app is free and included with 3Shape TRIOS and TRIOS Design Studio, Splint Studio, and Implant Studio.

Model Maker includes:

- Automatic hollowing - optimizes resin use*
- 'Drain holes' ensures a good printing experience*
- 'Bevels' added to base structure to ensure models can be removed easily from the build platform

Model Maker is available in Europe as of today. The US and world – February 28.

Read more about Model Maker – <https://www.3shape.com/en/software-updates/model-maker>

Dental Practice Management Software integrations

The 3Shape Unite platform now integrates with Dental Practice Management Software. Dental professionals can setup their practice management integrations via the Unite Store. The integration enables professionals to automate their patient management via Unite and avoid workflow manual errors like filling in forms twice, save steps, and even open Unite directly from the practice management system for faster, more efficient workflows.

3Shape Automate - AI powered design platform for crowns & nightguards

Since launch, 3Shape Automate has delivered over 400,000 designs at a 90% acceptance rate to dental professionals. The AI powered 24/7 self-service design service provides dental designs for crowns and nightguards. Geared to dental lab technicians needing to cope with rush periods and

overflow orders, 3Shape Automate is a completely automated service. Dental technicians simply provide their material, manufacturing, contacts, occlusion preferences, and Automate sends back a dental design in as little as five minutes - created entirely by AI technology. The technician can accept the design, edit it in 3Shape Dental System, or instead, decline the design at no cost.

*The algorithm in Model Maker creates a model with hollowing and drain-holes. In some instances, i.e., when users do quadrant scans instead of full arch scans, the algorithm will choose to create a solid model to ensure mechanical stability.

About 3Shape

3Shape is changing dentistry together with dental professionals across the world by developing innovations that provide superior dental care for patients. We partner with industry leaders to give dental professionals open choices for their patient care as well as supporting professionals' continued education. 3Shape's solution portfolio includes the multiple award-winning 3Shape TRIOS® intraoral scanner and dental lab scanners, design services, and market-leading scanning and design software solutions for both dental practices and labs.

Our 1900 employees serve customers in over 100 countries from 3Shape offices around the world. 3Shape's products and innovations continue to challenge traditional methods, enabling dental professionals to treat more patients more effectively. www.3shape.com

